


Product Datasheet

EIF3H Antibody (orb1240966)

Description

EIF3H Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EIF3H.
Target	EIF3H
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store EIF3H antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	EIF3H antibody can be used for detection of EIF3H by ELISA at 1:12500. EIF3H antibody can be used for detection of EIF3H by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	O15372
NCBI	NP_003747
Dilution Range	EIF3H antibody can be used for detection of EIF3H by ELISA at 1:12500. EIF3H antibody can be used for detection of EIF3H by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Mouse Brain at 2 ...



Antibody used in WB on Human Muscle at 0...



See IHC 1 Data and Customer Feedback for more information

Antibody used in IHC on Mouse brain at: ...