

Product Datasheet

LPPR2 Antibody (orb1241014)


Description

LPPR2 Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LPPR2.
Target	LPPR2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store LPPR2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	LPPR2 antibody can be used for detection of LPPR2 by ELISA at 1:12500. LPPR2 antibody can be used for detection of LPPR2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q96GM1
NCBI	NP_073574
Dilution Range	LPPR2 antibody can be used for detection of LPPR2 by ELISA at 1:12500. LPPR2 antibody can be used for detection of LPPR2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HeLa at 0.2...



Antibody used in WB on Hum. Fetal Brain ...