

Product Datasheet

EOGT Antibody (orb1241042)

Description

EOGT Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Brain at 0....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C3orf64.

Target

EOGT

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C3orf64 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C3orf64 antibody can be used for detection of C3orf64 by ELISA at 1:312500. C3orf64 antibody can be used for detection of C3orf64 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

52 kDa

Uniprot ID

[Q5NDL2](#)

NCBI

[NP_775925](#)

Dilution Range

C3orf64 antibody can be used for detection of C3orf64 by ELISA at 1:312500. C3orf64 antibody can be used for detection of C3orf64 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.