

Product Datasheet

ACCS Antibody (orb1241061)

Description

ACCS Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at 2....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PHACS.

Target

ACCS

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PHACS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PHACS antibody can be used for detection of PHACS by ELISA at 1:62500. PHACS antibody can be used for detection of PHACS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[Q96QU6](#)

NCBI

[EAW68073](#)

Dilution Range

PHACS antibody can be used for detection of PHACS by ELISA at 1:62500. PHACS antibody can be used for detection of PHACS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.