

Product Datasheet

GADD45B Antibody (orb1241130)

Description


GADD45B Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GADD45B.
Target	GADD45B
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GADD45B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GADD45B antibody can be used for detection of GADD45B by ELISA at 1:1562500. GADD45B antibody can be used for detection of GADD45B by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	O75293
NCBI	NP_056490
Dilution Range	GADD45B antibody can be used for detection of GADD45B by ELISA at 1:1562500. GADD45B antibody can be used for detection of GADD45B by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.







Antibody used in WB on Human HepG2 at 0....



Antibody used in IHC on Human Prostate C...



Antibody used in IHC on Human Prostate C...