

Product Datasheet

LDHB Antibody (orb1241148)

Description

LDHB Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human LDHB.

Target

LDHB

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store LDHB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

LDHB antibody can be used for detection of LDHB by ELISA at 1:62500. LDHB antibody can be used for detection of LDHB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

37 kDa

Uniprot ID

[P07195](#)

NCBI


[NP_002291](#)

Dilution Range


LDHB antibody can be used for detection of LDHB by ELISA at 1:62500. LDHB antibody can be used for detection of LDHB by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.




Antibody used in IHC on Human Placenta a...



See Immunoblot Data and Customer Feedback for more information

Antibody used in WB on Human pancreas at...



Green: LDHB
Blue: DAPI

See IHC 2 Data and Customer Feedback for more information

Antibody used in IHC on Human Lung at: 1...