

Product Datasheet

Antp Antibody (orb1241193)

Description	Antp Antibody	
Species/Host	Rabbit	
Reactivity	Drosophila	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of Drosophila ANTP.	Antibody used in WB on Drosophila at 0.2...
Target	Antp	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ANTP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	ANTP antibody can be used for detection of ANTP by ELISA at 1:62500. ANTP antibody can be used for detection of ANTP by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	33 kDa	
Uniprot ID	Q7KSY7	
NCBI	NP_996176	
Dilution Range	ANTP antibody can be used for detection of ANTP by ELISA at 1:62500. ANTP antibody can be used for detection of ANTP by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	