

Product Datasheet

cad Antibody (orb1241199)

Description

cad Antibody

Species/Host

Rabbit

Reactivity

Drosophila

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Drosophila
at 0.2...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of Drosophila CAD.

Target

cad

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CAD antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

CAD antibody can be used for detection of CAD by ELISA at 1:12500. CAD antibody can be used for detection of CAD by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

[P09085](#)

NCBI

[NP_599128](#)

Dilution Range

CAD antibody can be used for detection of CAD by ELISA at 1:12500. CAD antibody can be used for detection of CAD by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.