

Product Datasheet

GAS7 Antibody (orb1241204)

Description

GAS7 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human GAS7.

Target

GAS7

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GAS7 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GAS7 antibody can be used for detection of GAS7 by ELISA at 1:312500. GAS7 antibody can be used for detection of GAS7 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

54 kDa, 39 kDa, 48 kDa

Uniprot ID

[Q7Z571](#)

NCBI


[NP_958839](#)

Dilution Range

GAS7 antibody can be used for detection of GAS7 by ELISA at 1:312500. GAS7 antibody can be used for detection of GAS7 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 5....