

Product Datasheet

PSMC3IP Antibody (orb1241219)

Description

PSMC3IP Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSMC3IP.
Target	PSMC3IP
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PSMC3IP antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PSMC3IP antibody can be used for detection of PSMC3IP by ELISA at 1:62500. PSMC3IP antibody can be used for detection of PSMC3IP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	Q9P2W1
NCBI	NP_057640
Dilution Range	PSMC3IP antibody can be used for detection of PSMC3IP by ELISA at 1:62500. PSMC3IP antibody can be used for detection of PSMC3IP by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 1....