

Product Datasheet

PA2G4 Antibody (orb1241221)

Description

PA2G4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PA2G4.

Target

PA2G4

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PA2G4 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PA2G4 antibody can be used for detection of PA2G4 by ELISA at 1:1562500. PA2G4 antibody can be used for detection of PA2G4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

44 kDa

Uniprot ID

[Q9UQ80](#)

NCBI



[NP_006182](#)

Dilution Range

PA2G4 antibody can be used for detection of PA2G4 by ELISA at 1:1562500. PA2G4 antibody can be used for detection of PA2G4 by western blot at 1 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.


Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.

PA2G4



Antibody used in WB on Hum. Fetal Brain ...



Antibody used in WB on Human HepG2 cells...