

Product Datasheet

MAFK Antibody (orb1241226)

Description	MAFK Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human MAFK.	Antibody used in WB on Human Jurkat 2.5 ...
Target	MAFK	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store MAFK antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	MAFK antibody can be used for detection of MAFK by ELISA at 1:312500. MAFK antibody can be used for detection of MAFK by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	17 kDa	
Uniprot ID	O60675	
NCBI	NP_002351	
Dilution Range	MAFK antibody can be used for detection of MAFK by ELISA at 1:312500. MAFK antibody can be used for detection of MAFK by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	