

Product Datasheet

TSPAN15 Antibody (orb1241273)

Description

TSPAN15 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TSPAN15.

Target

TSPAN15

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store TSPAN15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

TSPAN15 antibody can be used for detection of TSPAN15 by ELISA at 1:312500. TSPAN15 antibody can be used for detection of TSPAN15 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

33 kDa

Uniprot ID

[O95858](#)

NCBI

[NP_036471](#)

Dilution Range

TSPAN15 antibody can be used for detection of TSPAN15 by ELISA at 1:312500. TSPAN15 antibody can be used for detection of TSPAN15 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....