



Product Datasheet

NGAL antibody (orb12413)

Description	Rabbit polyclonal antibody to NGAL
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Lipocalin-2/NGAL(NP_005555.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20KDa
Uniprot ID	P80188
Dilution Range	WB: 1:500-1:2000, IF/ICC: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A-549 cel...



Western blot analysis of extracts of A-5...