

Product Datasheet

COX15 Antibody (orb1241319)

Description

COX15 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human COX15.

Target

COX15

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store COX15 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

COX15 antibody can be used for detection of COX15 by ELISA at 1:62500. COX15 antibody can be used for detection of COX15 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

44 kDa, 46 kDa

Uniprot ID

[Q7KZN9](#)

NCBI



[NP_004367](#)

Dilution Range


COX15 antibody can be used for detection of COX15 by ELISA at 1:62500. COX15 antibody can be used for detection of COX15 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Muscle at ...



See Immunoblot 2 Data and Customer Feedback for more information

Antibody used in WB on Human fibroblast ...



Antibody used in WB on Human HepG2 cells...