

Product Datasheet

OSMR Antibody (orb1241324)

Description

OSMR Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HT1080 at
0...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human OSMR.

Target

OSMR

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store OSMR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

OSMR antibody can be used for detection of OSMR by ELISA at 1:12500. OSMR antibody can be used for detection of OSMR by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

110 kDa

Uniprot ID

[Q99650](#)

NCBI

[NP_003990](#)

Dilution Range

OSMR antibody can be used for detection of OSMR by ELISA at 1:12500. OSMR antibody can be used for detection of OSMR by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.