

Product Datasheet

ST3GAL5 Antibody (orb1241325)

Description

ST3GAL5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ST3GAL5.

Target

ST3GAL5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store ST3GAL5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

ST3GAL5 antibody can be used for detection of ST3GAL5 by ELISA at 1:1562500. ST3GAL5 antibody can be used for detection of ST3GAL5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

48 kDa, 46 kDa

Uniprot ID

[Q6NZX4](#)

NCBI


[NP_003887](#)

Dilution Range


ST3GAL5 antibody can be used for detection of ST3GAL5 by ELISA at 1:1562500. ST3GAL5 antibody can be used for detection of ST3GAL5 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.25...