


Product Datasheet

PSME3 Antibody (orb1241354)


Description

PSME3 Antibody


Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PSME3.
Target	PSME3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PSME3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PSME3 antibody can be used for detection of PSME3 by ELISA at 1:62500. PSME3 antibody can be used for detection of PSME3 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	31 kDa, 30 kDa
Uniprot ID	P61289
NCBI	NP_789839
Dilution Range	PSME3 antibody can be used for detection of PSME3 by ELISA at 1:62500. PSME3 antibody can be used for detection of PSME3 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Kidney at ...



Antibody used in IHC on Human Tonsil at ...



Antibody used in WB on Human HepG2 cells...