

Product Datasheet

NOLC1 Antibody (orb1241361)

Description

NOLC1 Antibody

Species/Host

Rabbit

Reactivity

Canine, Drosophila, Human, Mouse, Plant, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human NOLC1.

Target

NOLC1

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store NOLC1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

NOLC1 antibody can be used for detection of NOLC1 by ELISA at 1:62500. NOLC1 antibody can be used for detection of NOLC1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

73 kDa

Uniprot ID

[Q14978](#)

NCBI


[NP_004732](#)

Dilution Range


NOLC1 antibody can be used for detection of NOLC1 by ELISA at 1:62500. NOLC1 antibody can be used for detection of NOLC1 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in IHC on Zebrafish at 0.0...



Antibody used in WB on Human HepG2 cells...