

Product Datasheet

EIF2S3 Antibody (orb1241376)

Description

EIF2S3 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 5.0 ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EIF2S3.

Target

EIF2S3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store EIF2S3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

EIF2S3 antibody can be used for detection of EIF2S3 by ELISA at 1:312500. EIF2S3 antibody can be used for detection of EIF2S3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

51 kDa

Uniprot ID

[P41091](#)

NCBI

[NP_001406](#)

Dilution Range

EIF2S3 antibody can be used for detection of EIF2S3 by ELISA at 1:312500. EIF2S3 antibody can be used for detection of EIF2S3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.