



Product Datasheet

RRP1 Antibody (orb1241399)

Description	RRP1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human RRP1.	Antibody used in IHC on Human kidney.
Target	RRP1	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store RRP1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	Antibody used in WB on Human HepG2 at 0....
Note	For research use only	
Application notes	RRP1 antibody can be used for detection of RRP1 by ELISA at 1:62500. RRP1 antibody can be used for detection of RRP1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	53 kDa	
Uniprot ID	P56182	Antibody used in WB on Human, Mouse at 1...
NCBI	NP_003674	
Dilution Range	RRP1 antibody can be used for detection of RRP1 by ELISA at 1:62500. RRP1 antibody can be used for detection of RRP1 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	