

Product Datasheet

PFKL Antibody (orb1241402)

Description

PFKL Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PFKL.

Target

PFKL

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PFKL antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PFKL antibody can be used for detection of PFKL by ELISA at 1:312500. PFKL antibody can be used for detection of PFKL by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

64 kDa, 90 kDa, 85 kDa

Uniprot ID

[Q7L2M7](#)

NCBI

[AAH18295](#)

Dilution Range


PFKL antibody can be used for detection of PFKL by ELISA at 1:312500. PFKL antibody can be used for detection of PFKL by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in WB on Human HepG2 at 2....



Antibody used in WB on Human, Mouse at 1...