


Product Datasheet

PDE9A Antibody (orb1241405)


Description

PDE9A Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PDE9A.
Target	PDE9A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PDE9A antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PDE9A antibody can be used for detection of PDE9A by ELISA at 1:312500. PDE9A antibody can be used for detection of PDE9A by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	O76083
NCBI	NP_001001581
Dilution Range	PDE9A antibody can be used for detection of PDE9A by ELISA at 1:312500. PDE9A antibody can be used for detection of PDE9A by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Mouse Brain at 1:...



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1 ug...