


Product Datasheet

SSR2 Antibody (orb1241450)


Description

SSR2 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SSR2.
Target	SSR2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SSR2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SSR2 antibody can be used for detection of SSR2 by ELISA at 1:312500. SSR2 antibody can be used for detection of SSR2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	20 kDa, 20 kDa
Uniprot ID	P43308
NCBI	NP_003136
Dilution Range	SSR2 antibody can be used for detection of SSR2 by ELISA at 1:312500. SSR2 antibody can be used for detection of SSR2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Muscle.



Antibody used in WB on Human Jurkat 1.25...