

Product Datasheet

C1QB Antibody (orb1241495)

Description

C1QB Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C1QB.

Target

C1QB

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C1QB antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C1QB antibody can be used for detection of C1QB by ELISA at 1:1562500. C1QB antibody can be used for detection of C1QB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

27 kDa

Uniprot ID

[P02746](#)

NCBI


[NP_000482](#)

Dilution Range

C1QB antibody can be used for detection of C1QB by ELISA at 1:1562500. C1QB antibody can be used for detection of C1QB by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Muscle.



Antibody used in WB on Human Jurkat 0.5 ...