

Product Datasheet


TGFBI Antibody (orb1241502)

Description


TGFBI Antibody

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human TGFBI.
Target	TGFBI
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store TGFBI antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	TGFBI antibody can be used for detection of TGFBI by ELISA at 1:312500. TGFBI antibody can be used for detection of TGFBI by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	46 kDa, 75 kDa
Uniprot ID	Q15582
NCBI	AAH26352
Dilution Range	TGFBI antibody can be used for detection of TGFBI by ELISA at 1:312500. TGFBI antibody can be used for detection of TGFBI by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.


TGFBI



Antibody used in IF on Human cornea at 1...



Antibody used in WB on Human Jurkat 1.25...



Antibody used in WB on Jurkat 2.5 ug/ml ...