


Product Datasheet

C8B Antibody (orb1241508)


Description

C8B Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C8B.
Target	C8B
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store C8B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	C8B antibody can be used for detection of C8B by ELISA at 1:62500. C8B antibody can be used for detection of C8B by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P07358
NCBI	EAX06643
Dilution Range	C8B antibody can be used for detection of C8B by ELISA at 1:62500. C8B antibody can be used for detection of C8B by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human HepG2 at 0....



Antibody used in IHC on Human Liver.