

Product Datasheet

C8B Antibody (orb1241509)

Description

C8B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human C8B.

Target

C8B

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store C8B antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

C8B antibody can be used for detection of C8B by ELISA at 1:312500. C8B antibody can be used for detection of C8B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

61 kDa

Uniprot ID

[P07358](#)

NCBI


[NP_000057](#)

Dilution Range


C8B antibody can be used for detection of C8B by ELISA at 1:312500. C8B antibody can be used for detection of C8B by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Fetal liver...