

Product Datasheet

SLC26A5 Antibody (orb1241513)

Description

SLC26A5 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC26A5.

Target

SLC26A5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC26A5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC26A5 antibody can be used for detection of SLC26A5 by ELISA at 1:312500. SLC26A5 antibody can be used for detection of SLC26A5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

81 kDa, 75 kDa, 56 kDa

Uniprot ID

[P58743](#)

NCBI


[NP_945350](#)

Dilution Range

SLC26A5 antibody can be used for detection of SLC26A5 by ELISA at 1:312500. SLC26A5 antibody can be used for detection of SLC26A5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human HepG2 at 2....