

Product Datasheet

SLC39A12 Antibody (orb1241522)

Description	SLC39A12 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	Antibody used in WB on Human Jurkat 0.2-...
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC39A12.	
Target	SLC39A12	
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC39A12 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	SLC39A12 antibody can be used for detection of SLC39A12 by ELISA at 1:62500. SLC39A12 antibody can be used for detection of SLC39A12 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	73 kDa	
Uniprot ID	Q504Y0	
NCBI	NP_689938	
Dilution Range	SLC39A12 antibody can be used for detection of SLC39A12 by ELISA at 1:62500. SLC39A12 antibody can be used for detection of SLC39A12 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	