

Product Datasheet

SLC2A10 Antibody (orb1241533)

Description

SLC2A10 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
Jurkat 2.5 ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC2A10.

Target

SLC2A10

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC2A10 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC2A10 antibody can be used for detection of SLC2A10 by ELISA at 1:312500. SLC2A10 antibody can be used for detection of SLC2A10 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

60 kDa

Uniprot ID

[O95528](#)

NCBI

[NP_110404](#)

Dilution Range

SLC2A10 antibody can be used for detection of SLC2A10 by ELISA at 1:312500. SLC2A10 antibody can be used for detection of SLC2A10 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.