

Product Datasheet

SLC35F5 Antibody (orb1241534)

Description

SLC35F5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on
Transfected
293T ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC35F5.

Target

SLC35F5

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC35F5 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC35F5 antibody can be used for detection of SLC35F5 by ELISA at 1:1562500. SLC35F5 antibody can be used for detection of SLC35F5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

59 kDa

Uniprot ID

[Q8WV83](#)

NCBI

[NP_079457](#)

Dilution Range

SLC35F5 antibody can be used for detection of SLC35F5 by ELISA at 1:1562500. SLC35F5 antibody can be used for detection of SLC35F5 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.