


## Product Datasheet

### SLC2A6 Antibody (orb1241543)

## Description

SLC2A6 Antibody



### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, IHC, WB

### Applications

### Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC2A6.

### Target

SLC2A6

### Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC2A6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

### Note

For research use only

### Application notes

SLC2A6 antibody can be used for detection of SLC2A6 by ELISA at 1:312500. SLC2A6 antibody can be used for detection of SLC2A6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Clonality

Polyclonal

### MW

56 kDa

### Uniprot ID

[Q9UGQ3](#)

### NCBI


[NP\\_060055](#)

### Dilution Range


SLC2A6 antibody can be used for detection of SLC2A6 by ELISA at 1:312500. SLC2A6 antibody can be used for detection of SLC2A6 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

### Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human Muscle.



Antibody used in WB on Human Jurkat 2.5 ...