

Product Datasheet

SLC39A6 Antibody (orb1241551)

Description

SLC39A6 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC39A6.

Target

SLC39A6

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC39A6 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

SLC39A6 antibody can be used for detection of SLC39A6 by ELISA at 1:1562500. SLC39A6 antibody can be used for detection of SLC39A6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

48 kDa, 85 kDa

Uniprot ID

[Q13433](#)

NCBI


[NP_001092876](#)

Dilution Range


SLC39A6 antibody can be used for detection of SLC39A6 by ELISA at 1:1562500. SLC39A6 antibody can be used for detection of SLC39A6 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in WB on Human Jurkat 1.25...



Antibody used in WB on Hum. Fetal Brain ...