

Product Datasheet

SLC18A2 Antibody (orb1241560)

Description

SLC18A2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SLC18A2.

Antibody used in WB on Human Jurkat 0.2-...

Target

SLC18A2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store SLC18A2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on mouse and rat br...

Note

For research use only

Application notes

SLC18A2 antibody can be used for detection of SLC18A2 by ELISA at 1:312500. SLC18A2 antibody can be used for detection of SLC18A2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

56 kDa

Uniprot ID

[Q05940](#)

NCBI

[NP_003045](#)

Dilution Range

SLC18A2 antibody can be used for detection of SLC18A2 by ELISA at 1:312500. SLC18A2 antibody can be used for detection of SLC18A2 by western blot at 0.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.