

Product Datasheet

STEAP3 Antibody (orb1241591)

Description

STEAP3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Antibody
used in WB
on Human
HepG2 at
5....

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human STEAP3.

Target

STEAP3

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store STEAP3 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

STEAP3 antibody can be used for detection of STEAP3 by ELISA at 1:312500. STEAP3 antibody can be used for detection of STEAP3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

54 kDa, 55 kDa, 56 kDa

Uniprot ID

[Q658P3](#)

NCBI

[NP_001008410](#)

Dilution Range

STEAP3 antibody can be used for detection of STEAP3 by ELISA at 1:312500. STEAP3 antibody can be used for detection of STEAP3 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.