

Product Datasheet


M6PR Antibody (orb1241594)

Description

M6PR Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human M6PR.
Target	M6PR
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store M6PR antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	M6PR antibody can be used for detection of M6PR by ELISA at 1:312500. M6PR antibody can be used for detection of M6PR by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	P20645
NCBI	NP_002346
Dilution Range	M6PR antibody can be used for detection of M6PR by ELISA at 1:312500. M6PR antibody can be used for detection of M6PR by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.

M6PR



See Unconjugated 20kDa and Customer Feedback for more information

Antibody used in WB on Human, Mouse cell...



Antibody used in WB on Human Jurkat 5.0 ...