

Product Datasheet

PRPF19 Antibody (orb1241629)

Description

PRPF19 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PRPF19.

Target

PRPF19

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PRPF19 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PRPF19 antibody can be used for detection of PRPF19 by ELISA at 1:1562500. PRPF19 antibody can be used for detection of PRPF19 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

55 kDa

Uniprot ID

[Q9UMS4](#)

NCBI

[NP_055317](#)


Dilution Range

PRPF19 antibody can be used for detection of PRPF19 by ELISA at 1:1562500. PRPF19 antibody can be used for detection of PRPF19 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

Antibody used in WB on Human Jurkat 2.5 ...



Antibody used in WB on mouse brain at 1:...