

Product Datasheet

PANX1 Antibody (orb1241670)


Description

PANX1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PANX1.
Target	PANX1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PANX1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PANX1 antibody can be used for detection of PANX1 by ELISA at 1:62500. PANX1 antibody can be used for detection of PANX1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q96RD7
NCBI	NP_056183
Dilution Range	PANX1 antibody can be used for detection of PANX1 by ELISA at 1:62500. PANX1 antibody can be used for detection of PANX1 by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in IHC on mouse retina and...



Antibody used in WB on Human Jurkat 0.2-...