

Product Datasheet

PANX2 Antibody (orb1241671)

Description

PANX2 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PANX2.

Target

PANX2

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store PANX2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

PANX2 antibody can be used for detection of PANX2 by ELISA at 1:62500. PANX2 antibody can be used for detection of PANX2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

70 kDa

Uniprot ID

[Q96RD6](#)

NCBI


[NP_443071](#)

Dilution Range

PANX2 antibody can be used for detection of PANX2 by ELISA at 1:62500. PANX2 antibody can be used for detection of PANX2 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date


12 months from date of receipt.



Antibody used in IHC on Human Muscle.



Antibody used in WB on Human HepG2 at 1....



Antibody used in WB on HepG2 at 2.5 µg/mL...