

Product Datasheet

GNAS Antibody (orb1241679)

Description

GNAS Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GNAS.

Target

GNAS

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store GNAS antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Note

For research use only

Application notes

GNAS antibody can be used for detection of GNAS by ELISA at 1:312500. GNAS antibody can be used for detection of GNAS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

42 kDa

Uniprot ID

[Q5FWY2](#)

NCBI


[NP_536351](#)

Dilution Range


GNAS antibody can be used for detection of GNAS by ELISA at 1:312500. GNAS antibody can be used for detection of GNAS by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.



Antibody used in IHC on Human kidney.



Antibody used in IHC on Human Lung.



Antibody used in WB on Human Jurkat 2.5 ...