


Product Datasheet

S1PR5 Antibody (orb1241685)


Description

S1PR5 Antibody

Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human EDG8.
Target	S1PR5
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store EDG8 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	EDG8 antibody can be used for detection of EDG8 by ELISA at 1:312500. EDG8 antibody can be used for detection of EDG8 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	16 kDa, 17 kDa, 14 kDa
Uniprot ID	Q9H228
NCBI	EAW84115
Dilution Range	EDG8 antibody can be used for detection of EDG8 by ELISA at 1:312500. EDG8 antibody can be used for detection of EDG8 by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Intestine.



Antibody used in WB on Human Jurkat 2.5 ...