



Product Datasheet

MMP9 antibody (orb12417)


Description	Rabbit polyclonal antibody to MMP9
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human MMP9(NP_004985.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	92KDa
Uniprot ID	P14780
Dilution Range	WB: 1:200-1:500, IHC-P: 1:100-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunohistochemical staining of human es...



Immunohistochemical staining of human li...