


Product Datasheet

BPIFA1 Antibody (orb1241712)


Description

BPIFA1 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human PLUNC.
Target	BPIFA1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store PLUNC antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	PLUNC antibody can be used for detection of PLUNC by ELISA at 1:1562500. PLUNC antibody can be used for detection of PLUNC by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	28 kDa, 27 kDa
Uniprot ID	Q9NP55
NCBI	NP_057667
Dilution Range	PLUNC antibody can be used for detection of PLUNC by ELISA at 1:1562500. PLUNC antibody can be used for detection of PLUNC by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in WB on Human Jurkat 0.2-...



Antibody used in IHC on Human Lung at 5....