


Product Datasheet

GAPVD1 Antibody (orb1241718)


Description

GAPVD1 Antibody

Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human GAPVD1.
Target	GAPVD1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store GAPVD1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	GAPVD1 antibody can be used for detection of GAPVD1 by ELISA at 1:312500. GAPVD1 antibody can be used for detection of GAPVD1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	166 kDa
Uniprot ID	Q14C86
NCBI	NP_056450
Dilution Range	GAPVD1 antibody can be used for detection of GAPVD1 by ELISA at 1:312500. GAPVD1 antibody can be used for detection of GAPVD1 by western blot at 0.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Brain.



Antibody used in WB on Human HepG2 at 0....