


Product Datasheet

CTDSPL Antibody (orb1241742)

Description

CTDSPL Antibody



Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody used in IHC on Human Kidney at ...

Immunogen

Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human CTDSPL.

Target

CTDSPL

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store CTDSPL antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in WB on Human DLD1 cells ...

Note

For research use only

Application notes

CTDSPL antibody can be used for detection of CTDSPL by ELISA at 1:312500. CTDSPL antibody can be used for detection of CTDSPL by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

30 kDa

Uniprot ID

[Q3ZTU0](#)

NCBI

[NP_005799](#)

Dilution Range

CTDSPL antibody can be used for detection of CTDSPL by ELISA at 1:312500. CTDSPL antibody can be used for detection of CTDSPL by western blot at 1.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.