


Product Datasheet

SEMG1 Antibody (orb1241774)


Description

SEMG1 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human SEMG1.
Target	SEMG1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store SEMG1 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	SEMG1 antibody can be used for detection of SEMG1 by ELISA at 1:62500. SEMG1 antibody can be used for detection of SEMG1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	51 kDa, 45 kDa
Uniprot ID	P04279
NCBI	NP_002998
Dilution Range	SEMG1 antibody can be used for detection of SEMG1 by ELISA at 1:62500. SEMG1 antibody can be used for detection of SEMG1 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.



Antibody used in IHC on Human Lung at 4....



Antibody used in IHC on Human Kidney at ...



Antibody used in WB on Human Jurkat 1.25...