

Product Datasheet

ALAS2 Antibody (orb1241805)

Description


ALAS2 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human ALAS2.
Target	ALAS2
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store ALAS2 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.
Note	For research use only
Application notes	ALAS2 antibody can be used for detection of ALAS2 by ELISA at 1:62500. ALAS2 antibody can be used for detection of ALAS2 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Clonality	Polyclonal
MW	57 kDa, 65 kDa, 61 kDa, 63 kDa
Uniprot ID	P22557
NCBI	NP_001033058
Dilution Range	ALAS2 antibody can be used for detection of ALAS2 by ELISA at 1:62500. ALAS2 antibody can be used for detection of ALAS2 by western blot at 5.0 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.
Expiration Date	12 months from date of receipt.






Antibody used in WB on Human Jurkat 5.0 ...



Antibody used in IHC on Human Kidney at ...