

Product Datasheet

STK11 Antibody (orb1241831)

Description

STK11 Antibody

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat, Zebrafish

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Antibody produced in rabbits immunized with a synthetic peptide corresponding to a region of human STK11.

Antibody used in WB on Human Jurkat 1.25...

Target

STK11

Preservatives

Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

For short periods of storage (days) store at 4°C. For longer periods of storage, store STK11 antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.

Antibody used in IHC on Human Skeletal m...

Note

For research use only

Application notes

STK11 antibody can be used for detection of STK11 by ELISA at 1:1562500. STK11 antibody can be used for detection of STK11 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Clonality

Polyclonal

MW

48 kDa

Uniprot ID

[Q15831](#)

NCBI

[NP_000446](#)

Dilution Range

STK11 antibody can be used for detection of STK11 by ELISA at 1:1562500. STK11 antibody can be used for detection of STK11 by western blot at 1.25 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.

Expiration Date

12 months from date of receipt.

