


Product Datasheet

UROD Antibody (orb1241836)

Description	UROD Antibody	
Species/Host	Rabbit	 Antibody used in IHC on Human Kidney at ...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Antibody produced in rabbits immunized with a synthetic peptide corresponding a region of human UROD.	
Target	UROD	 Antibody used in WB on Human HepG2 cells...
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	For short periods of storage (days) store at 4°C. For longer periods of storage, store UROD antibody at -20°C. As with any antibody avoid repeat freeze-thaw cycles.	
Note	For research use only	
Application notes	UROD antibody can be used for detection of UROD by ELISA at 1:62500. UROD antibody can be used for detection of UROD by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Clonality	Polyclonal	
MW	40 kDa	
Uniprot ID	P06132	
NCBI	NP_000365	
Dilution Range	UROD antibody can be used for detection of UROD by ELISA at 1:62500. UROD antibody can be used for detection of UROD by western blot at 2.5 µg/mL, and HRP conjugated secondary antibody should be diluted 1:50,000 - 100,000.	
Expiration Date	12 months from date of receipt.	